

# Clean Air Filter

The Clean Air Filter effectively eliminates excess halogenated anesthetic gases after patient inhalation in the operating-room environment. Each canister offers an average useful life of 12 to 15 hours. It can remove no fewer than 50 grams of waste gases and then can be conveniently discarded. Not for use with flammable gases.

ITEM #	DESCRIPTION	UNIT
617029	Clean Air Filter	Each



Available exclusively through your veterinarian from:



(800) 824-3703  
www.mwivet.com



(888) 722-2242  
www.aahamarketlink.com