





Collar Instructions

VetOne® collars are designed to provide a comfortable, yet effective fit on both dogs and cats. They provide vital protection behind the neck of the patient so that wounds can rapidly heal without the patient licking or biting the wound area.

To apply the VetOne collar, please refer to the instruction chart below. To further assist you in selecting the correct size for the dog or cat, refer to the sizing chart. NOTE: Collar depth is measured from the outside edge to the inside edge.



SIZE	NECK CIRCUMFERENCE	COLLAR DEPTH
7.5	8-3/4" – 10-1/8"	3"
10	9" – 10-1/4"	4"
12	11" – 12-1/4"	5"
15	13" – 16"	6"
20	13-1/8" – 16"	8"
25	15-1/2" – 19"	9-7/8"
30	18-1/4" – 22-3/4"	11-7/8"