

# Surgical Adhesive

VetOne Surgical Adhesive contains purified n-butyl cyanoacrylate with a D&C-certified dye that allows for differentiation between the product and bodily fluid, but it dries clear for a clean appearance. It sets up quickly (5-15 seconds) to form a thin, flexible waterproof bandage that is sloughed naturally in 5–7 days. Disposable applicator tips (included) fit snugly into the tip on the bottle and allow greater control over the drop size and placement of the product.

VetOne Surgical Adhesive is designed for wound management in situations such as:

- Lacerations and abrasions
- Cat declaws
- Cuts
- Tail docking

For topical application only. This product is not intended for internal closures or for use on infected or deep puncture wounds.

ITEM#	DESCRIPTION	SIZE
510077	Surgical Adhesive	2 mL



Available exclusively through your veterinarian from:



(800) 824-3703  
www.mwivet.com



(888) 722-2242  
www.aahamarketlink.com