

Porta-Spin™ Portable Centrifuge

The Porta-Spin is designed to provide years of dependable service either in the field or in the lab; use it wherever a 12V DC outlet can be found, or take it back to use in the lab with the 120V AC adapter (sold separately). The fixed-angle 6-place rotor is designed for use with standard 2–15 mL tubes. It will also spin up to three 15 mL tubes. You can spin 2 mL and 4-5 mL tubes with adapters available from MWI. With 6 tubes in place, the nominal speed of the Porta-Spin is 3,700 RPM with the 12V DC power source and 3,500 RPM with the 120V AC power source. The side-mounted Tube Holdster™ rack (sold separately) keeps up to 18 tubes close at hand.

- Dependable mechanical 30-minute timer with bell signal at end of the cycle
- Fish-Eye Level built into the lid to make sure the unit is level before operation
- Safety switch in lid shuts off motor when lid is unlatched
- Latch release is located on the right front side of unit
- Indicator light illuminates when unit is powered on
- Sturdy stainless steel handle for easy carrying

ITEM#	DESCRIPTION	UNIT
619015	Port-Spin Centrifuge	Each
619016	Tube Holdster	Each
619017	120V AC Power Adapter	Each
605015	Adapters for 2 mL or Smaller Tubes	Pack of 2
606015	Adapters for 3.5 mL–7 mL Tubes	Pack of 10



Shown with optional Tube Holdster™.

Available exclusively through your veterinarian from:



(800) 824-3703
www.mwivet.com



(888) 722-2242
www.aahamarketlink.com